



Recombinant Salmonella typhimurium Enoyl-[acyl-carrier-protein] reductase [NADH] FabI (fabI)

| | |
|--------------------------|---|
| Product Code | CSB-EP325843SXB |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P16657 |
| Product Type | Recombinant Protein |
| Immunogen Species | Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | GFLSGKRIL VTGVASKLSI AYGIAQAMHR EGAELAFTYQ NDKLKGRVEE FAAQLGSSIV LPCDVAEDAS IDAMFAELGN VWPKFDGDFVH SIGFAPGDQL DGDYVNAVTR EGFKVAHDIS SYSFVAMAKA CRTMLNPGSA LLTLSYLGAE RAIPNYNVMG LAKASLEANV RYMANAMGPE GVRVNAISAG PIRTLAASGI KDFRKMLAHC EAVTPIRRTV TIEDVGNSAA FLCSDL SAGI SGEVVHVDGG FSIAAMNELE LK |
| Source | E.coli |
| Target Names | fabI |
| Protein Names | Recommended name: Enoyl-[acyl-carrier-protein] reductase [NADH] FabI Short name= ENR EC= 1.3.1.9 Alternative name(s): NADH-dependent enoyl-ACP reductase |
| Expression Region | 2-262 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |