



Recombinant Salmonella typhimurium Enoyl-[acyl-carrier-protein] reductase [NADH] FabI (fabI)

Product Code	CSB-MP325843SXB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16657
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	≥85% (SDS-PAGE)
Sequence	GFLSGKRIL VTGVASKLSI AYGIAQAMHR EGAELAFTYQ NDKLKGRVEE FAAQLGSSIV LPCDVAEDAS IDAMFAELGN VWPKFDGFMVH SIGFAPGDQL DGDYVNAVTR EGFKVAHDIS SYSFVAMAKA CRTMLNPGSA LLTLSYLGAE RAIPNYNVMG LAKASLEANV RYMANAMGPE GVRVNAISAG PIRTLAASGI KDFRKMLAHC EAVTPIRRTV TIEDVGNSAA FLCSDL SAGI SGEVVHVDGG FSIAAMNELE LK
Source	Mammalian cell
Target Names	fabI
Protein Names	Recommended name: Enoyl-[acyl-carrier-protein] reductase [NADH] FabI Short name= ENR EC= 1.3.1.9 Alternative name(s): NADH-dependent enoyl-ACP reductase
Expression Region	2-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.