



Recombinant *Shigella flexneri* Invasin IpaD (ipaD)

Product Code	CSB-YP325954SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18013
Product Type	Recombinant Protein
Immunogen Species	<i>Shigella flexneri</i>
Purity	≥85% (SDS-PAGE)
Sequence	MNITTLTNSI STSSFSPNNT NGSSTETVNS DIKTTTSSHP VSSLTMLNDT LHNIRTTNQA LKKELSQKTL TKTSLEEIAL HSSQISMDVN KSAQLLDILS RNEYPIKDA RELLHSAPKE AELDGDQMIS HRELWAKIAN SINDINEQYL KVYEHAVSSY TQMYQDFSAV LSSLAGWISP GGNDGNSVKL QVNSLKKALE ELKEKYKDKP LYPANNTVSQ EQANKWLTEL GGTIGKVSQK NGGYVVSINM TPIDNMLKSL DNLGGNGEVV LDNAKYQAWN AGFSAEDEM KNNLQTLVQK YSNANSIFDN LVKVLSSIS SCTDTDKLFL HF
Source	Yeast
Target Names	ipaD
Protein Names	Recommended name: Invasin ipaD Alternative name(s): 36 kDa membrane antigen
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.