



# Recombinant *Drosophila melanogaster* Protein tailless (tll)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP325960DLU-B   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | P18102  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Drosophila melanogaster</i> (Fruit fly)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MQSSEGSPDM MDQKYNVRL SPAASSRILY HVPCKVCRDH SSGKHYGIYA<br>CDGCAGFFKR SIRRSRQYVC KSQKQGLCVV DKTHRNQCRA CRLRKCFEVG<br>MNKDAVQHER GPRNSTLRRH MAMYKDAMMG AGEMPQIPAE<br>ILMNTAALTG FPGVPMMPG LPQRAGHHPA HMAAFQPPPS AAVALDLSVP<br>RVPHHPVHQQ HHGFFSPTAA YMNALATRAL PPTPLMAAE HIKETAAEHL<br>FKNVNWIKSV RAFTELPMPD QLLLLLEESWK EFFILAMAQY LMPMNFAQLL<br>FVYESENANR EIMGMTREV HAFQEVNLQL CHLNIDSTEY ECLRAISLFR<br>KSPPSASSTE DLANSSILTG SGSPNSSASA ESRGLLESGK VAAMHNDARS<br>ALHNYIQRTH PSQPMRFQTL LGVVQLMHKV SSFTIEELFF RKTIGDITIV<br>RLISDMYSQR KI |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | tll   |
| <b>Protein Names</b>     | Recommended name: Protein tailless Alternative name(s): Nuclear receptor subfamily 2 group E member 2   |
| <b>Expression Region</b> | 1-452   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |