



Recombinant Lysinibacillus sphaericus 8-amino-7-oxononanoate synthase (bioF)

Product Code	CSB-EP326014LRU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22806
Product Type	Recombinant Protein
Immunogen Species	Lysinibacillus sphaericus (Bacillus sphaericus)
Purity	≥85% (SDS-PAGE)
Sequence	MNDRFRRELQ VIEEQGLTRK LRLFSTGNES EVVMNGKKFL LFSSNNYLGL ATDSRLKKKA TEGISKYGTG AGGSRLTTGN FDIHEQLESE IADFKKTEAA IVFSSGYLAN VGVISSVMKA GDTIFSDAWN HASIIDGCRL SKAKTIVYEH ADMVDLERKL RQSHGDGLKF IVTDGVFSMD GDIAPLPKIV ELAKEYKAYI MIDDAHATGV LGNDGCGTAD YFGLKDEIDF TVGTLSKAIG AEGGFVSTSS IAKNYLLNNA RSFIFQTALS PSAIEAAREG ISIIQNEPER RKQLLKNAQY LRLKLEESGF VMKEGETPII SLIIGGSHEA MQFSAKLLDE GVFIPAIRPP TVPKGSSRLR ITVMATHTIE QLDMVISIKI KIGKEMGIV
Source	E.coli
Target Names	bioF
Protein Names	Recommended name: 8-amino-7-oxononanoate synthase Short name= AONS EC= 2.3.1.47 Alternative name(s): 7-keto-8-amino-pelargonic acid synthase Short name= 7-KAP synthase 8-amino-7-ketopelargonate synthase L-alanine-- pime
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.