



# Recombinant *Saccharomyces cerevisiae* Protein MAK32 (MAK32)

<b>Product Code</b>	CSB-MP326032SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P23060
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMNEEDSTET KSLVITNGMF IIDDIERSKY NIHYKNVPGG GGTFAILGAC IISGNVTSK GLKWIVDRGS DFPKEVIREI DSWGTDVRFRR DDFSRLTTKG LNYEYEGSDDL RKFKFLTPKK QINVDDWIST FGQKIIDEMH AFHLLCSGSR CLDIINDLLR VKSSKGTKPI VIWEPFPDLC DFDHQNDIKS VMQRNDVTVI LSPNAEESSR LFGLSSKEPT SLEECLALAH RFDDFMDENN MCILRCGALG SISVSEKFKN GRTYDHFYAY HFKTQSKVLD PTGGGNSFLG GFAVSIALTK SLDIASICGN IAAGAIIEQF GIPRYDPIAK TWNGITFLDR LKFYLSQSGL QYNINDLYKS LTR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	MAK32
<b>Protein Names</b>	Recommended name: Protein MAK32 Alternative name(s): Maintenance of killer protein 32
<b>Expression Region</b>	1-363
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.