



Recombinant *Rhodobacter capsulatus* 1-phosphofructokinase (fruK)

Product Code	CSB-BP326058RLC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23386
Product Type	Recombinant Protein
Immunogen Species	<i>Rhodobacter capsulatus</i> (<i>Rhodopseudomonas capsulata</i>)
Purity	>85% (SDS-PAGE)
Sequence	MTLRIATVSL NSAVDQTVTV PGFTADAVNR VAASRIDAGG KGVNVASFLA HVGHGVAVTG LLGAENAALF ARHFAATGLV DACQRLPGAT RTNVKIVDPL QDQVTDLNFP GIAAGPADLD AVAATLTELL AQGLDWVALC GSLPAGIGAE AYAELAALAR KGGARVALDT SGPALGLALA ARPDIVKPNV AELGAHLGRT LTGLESVREA ARDLAASGVG LVAVSMGAGG AVLVRGAEAV LAIPPATPIA STVGAGDAMV AGLIHAATLG LDLAETARLA TAFSLGALGE IGPLPPPER LAALARTVTV KTLPPV
Source	Baculovirus
Target Names	fruK
Protein Names	Recommended name: 1-phosphofructokinase EC= 2.7.1.56 Alternative name(s): Fructose 1-phosphate kinase
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.