



Recombinant *Saccharomyces cerevisiae* Proteasome component Y13 (PRE9)

Product Code	CSB-MP326080SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23638
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGSRRYDSRT TIFSPEGRLY QVEYALE SIS HAGTAIGIMA SDGIVLAAER KVTSTLLEQD TSTEKLYKLN DKIAVAVAGL TADAEILINT ARIHAQNYLK TYNEDIPVEI LVRRLSDIKQ GYTQHGGLRP FGVSFIYAGY DDRYGYQLYT SNPSGNYTGW KAISVGANTS AAQTLQMDY KDDMKVDDAI ELALKTLSKT TDSSALTYDR LEFATIRKGA NDGEVYQKIF KPQEIKDILV KTGITKKDED EEADEDMDK
Source	Mammalian cell
Target Names	PRE9
Protein Names	Recommended name: Proteasome component Y13 EC= 3.4.25.1 Alternative name(s): Macropain subunit Y13 Multicatalytic endopeptidase complex subunit Y13 Proteinase YSCE subunit 13
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.