



Recombinant Escherichia coli Transcriptional regulatory protein PhoP (phoP)

Product Code	CSB-YP326097ENV
Relevance	Member of the two-component regulatory system PhoP/PhoQ involved in adaptation to low Mg ²⁺ environments and the control of acid resistance genes. In low periplasmic Mg ²⁺ , PhoQ phosphorylates PhoP, resulting in the expression of PhoP-activated genes (PAG) and repression of PhoP-repressed genes (PRG). In high periplasmic Mg ²⁺ , PhoQ dephosphorylates phospho-PhoP, resulting in the repression of PAG and may lead to expression of some PRG. Mediates magnesium influx to the cytosol by activation of MgtA. Promotes expression of the two-component regulatory system rstA/rstB and transcription of the hemL, mgrB, nagA, slyB, vboR and yrbL genes.
Abbreviation	Recombinant E.coli phoP protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23836
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MRVLVVEDNALLRHHLKVQIQDAGHQVDDAEDAKEYDYLLNEHIPDIAIVDLGL PDEDGLSLIRRWRSNDVSLPILVLTARES WQDKVEVLSAGADDYVTKPFHIEEV MARMQALMRRNSGLASQVISLPPFQVDLSRRELSINDEVIKLTAFEYTIMETLIR NNGKVVS KDSLMLQLYPDAELRESHTIDVLMGRLRKKIQAQYPQEVITTVRGQ GYLFELR
Research Area	others
Source	Yeast
Target Names	phoP
Protein Names	Recommended name: Transcriptional regulatory protein phoP
Expression Region	1-223aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged and C-terminal Myc-tagged
Mol. Weight	29.0kDa
Protein Length	Full Length
Image	

