



Recombinant Escherichia coli Hydrogenase isoenzymes formation protein hypD (hypD)

Product Code	CSB-EP326123ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24192
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MRFVDEYRAP EQVMQLIEHL RERASHLSYT AERPLRIMEV CGGHTHAIFK FGLDQLLPEN VEFIHGPGCP VCVLPMGRID TCVEIASHPE VIFCTFGDAM RVPGKQGSLL QAKARGADVR IVYSPMDALK LAQENPTRKV VFFGLGFETT MPTTAITLQQ AKARDVQNFY FFCQHITLIP TLRSLLEQPD NGIDAFLAPG HVSMVIGTDA YNFIASDFHR PLVVAGFEPL DLLQGVVMLV QKIAAHSKV ENQYRRVVPD AGNLLAQQAI ADVFCVNGDS EWRGLGVIES SGVHLTPDYQ RFDAEAHFRP APQQVCDDPR ARCGEVLTGK CKPHQCPLFG NTCNPQTAFG ALMVSSEGAC AAWYQYRQQE SEA
Source	E.coli
Target Names	hypD
Protein Names	Recommended name: Hydrogenase isoenzymes formation protein hypD
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.