



# Recombinant Herpetosiphon aurantiacus Type-2 restriction enzyme HgiCI (hgiCIR)

<b>Product Code</b>	CSB-MP326195HGK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P25258
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Herpetosiphon aurantiacus (Herpetosiphon giganteus)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNYQRSFEDL EFNAIKWWPQ ELSATVAEAS VLPILISSQD LFISILKLSG THPEQIFDVI NAAQISANLF LKHLVVLADY GGEMIKRLGK EFQEIFPRMD STLEYMNYT FKGEQYTYIF KKLPIKGLDN SKLAIDGKAI IEIKPLSDLY RDMIMILLYG STTEQFNLAG LEKCEIGTIL GKNEIITYYI TQKYLYVSRI TNGANTNSLG QIAQTYVCDI LSKYLPNDYS VTRNGKILLS DLNSQDSTKT SFDILVEFAD KKVGEVFSFQ VTTNSTIERK AGQARDRQNR MAHAHYWIAW VIDGAGNFER SGAVRAICRY SDCTVAYSES EIAVLAAFIQ EKFN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	hgiCIR
<b>Protein Names</b>	Recommended name: Type-2 restriction enzyme HgiCI Short name= R.HgiCI EC= 3.1.21.4 Alternative name(s): Endonuclease HgiCI Type II restriction enzyme HgiCI
<b>Expression Region</b>	1-345
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.