



Recombinant Herpetosiphon aurantiacus Type-2 restriction enzyme HgiCI (hgiCIR)

Product Code	CSB-MP326195HGK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25258
Product Type	Recombinant Protein
Immunogen Species	Herpetosiphon aurantiacus (Herpetosiphon giganteus)
Purity	≥85% (SDS-PAGE)
Sequence	MNYQRSFEDL EFNAIKWWPQ ELSATVAEAS VLPILISSQD LFISILKLSG THPEQIFDVI NAAQISANLF LKHLVVLADY GGEMIKRLGK EFQEIFPRMD STLEYYMNYT FKGEQYTYIF KKLPIKGLDN SKLAIDGKAI IEIKPLSDLY RDMIMILLYG STTEQFNLAG LEKCEIGTIL GKNEIITYYI TQKYLYVSRI TNGANTNSLG QIAQTYVCDI LSKYLPNDYS VTRNGKILLS DLNSQDSTKT SFDILVEFAD KKVGEVFSFQ VTTNSTIERK AGQARDRQNR MAHAHYWIAW VIDGAGNFER SGAVRAICRY SDCTVAYSES EIAVLAAFIQ EKFN
Source	Mammalian cell
Target Names	hgiCIR
Protein Names	Recommended name: Type-2 restriction enzyme HgiCI Short name= R.HgiCI EC= 3.1.21.4 Alternative name(s): Endonuclease HgiCI Type II restriction enzyme HgiCI
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.