



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YCR102C (YCR102C)

Product Code	CSB-EP326229SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25608
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKAVVIEDGK AVVKEGVPIP ELEEGFVLIK TLAVAGNPTD WAHIDYKVG QGSILGCDAA GQIVKLGPAV DPKDFSIGDY IYGFIHGSSV RFPSNGAF YSAISTVVAY KSPNELKFLG EDVLPAGPVR SLEGAATIPV SLTTAGLV YNLGLNLKWE PSTPQRNGPI LLWGGATAVG QSLIQLANKL NGFTKIIV SRKHEKLLKE YGADQLFDYH DIDVVEQIKH KYNNISYLV CVANQNTLQ VYKCAADKQD ATVVELTNLT EENVKKNRR QNVTIDRTRL YSIGGHEVP GGITFPADPE ARRAATEFVK FINPKISDGQ IHHIPARVYK NGLYDVPRIL EDIKIGKNSG EKLVAVLN
Source	E.coli
Target Names	YCR102C
Protein Names	Recommended name: Uncharacterized protein YCR102C
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.