



Recombinant Arabidopsis thaliana Glyceraldehyde-3-phosphate dehydrogenase B, chloroplastic (GAPB)

Product Code	CSB-MP326249DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25857
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	KLKVAINGFG RIGRNFLRCW HGRKDSPLEV VVLNDSGGVK NASHLLKYDS MLGTFKAEVK IVDNETISVD GKLIKVVSNR DPLKLPWAEL GIDIVIEGTG VFVDGPGAGK HIQAGASKVI ITAPAKGADI PTYVMGVNEQ DYGHDEVANII SNASCTTNCL APFAKVLDEE FGIVKGTMTT THSYTGDQRL LDASHRDLRR ARAAALNIVP TSTGAAKAVS LVLPQLKGKL NGIALRVPTP NVSVVDLVIN VEKKGLTAED VNEAFRKAAN GPMKGILDVC DAPLVSVDFR CSDVSTTIDS SLTMVMGDDM VKVVAWYDNE WGYSQRVVDL AHLVASKWPG AEAVGSGDPL EDFCKTNPAD EECKVYD
Source	Mammalian cell
Target Names	GAPB
Protein Names	Recommended name: Glyceraldehyde-3-phosphate dehydrogenase B, chloroplastic EC= 1.2.1.13 Alternative name(s): NADP-dependent glyceraldehydephosphate dehydrogenase subunit B
Expression Region	81-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.