



# Recombinant 6-phospho-beta-glucosidase (arbB)

<b>Product Code</b>	CSB-MP326280EGI
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P26206
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Dickeya chrysanthemi (Pectobacterium chrysanthemi) (Erwinia chrysanthemi)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSNPFFPAHFL WGGAIANQV EGAYLTDGKG LSTSDLQPQG IFGEIVTRQP GDSGIKDVAI DFYHRYPQDI ALFAEMGFTC LRISIAWTRI FPQGDEAEPN EAGLAFYDRL FDELAKYGIQ PLVTLSHYEM PYGLVEKHGG WGNRLTIDCF ERYARTVFAR YRHKVKRWLT FNEINMSLHA PFTGVGLPPD SDKAAIQAI HHQLVASARA VKACHDMIPD AQIGNMLLGA MLYPLTSKPE DVMESLHQNR EWLFFGDVQV RGAYPGYMHR YFREQGITLN ITAQDKQDLK ATVDFISFSY YMTGCVTTDE AQLEKTRGNI LNMVNPYLE SSEWGWQIDP LGLRYLLNFL YDRYQKPLFI VENGLGAKDK IEENGDIYDD YRIRYLNHDL VQVGEAIDDG VEVLGYTCWG PIDLVSASKA EMSKRYGFIY VDRDDAGHGS LERRRKKSFY WYQSVIASHG KTLTR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	arbB
<b>Protein Names</b>	Recommended name: 6-phospho-beta-glucosidase EC= 3.2.1.86
<b>Expression Region</b>	1-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.