



Recombinant Pig Myocyte-specific enhancer factor 2C (MEF2C)

Product Code	CSB-MP013673PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4UTP7
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	MGRKKIQITR IMDERNRQVT FTKRKFGMLK KAYELSVLCD CEIALIIFNS TNKLFQYAST DMDKVLLKYT EYNEPHESRT NSDIVEALNK KENKGCESPD PDSSYALTPR TEEKYKKINE EFDNMIKSHK IPAVPPPNFE MPVSIPVSSH NSLVYSNPVS SLGNPNFLPL AHPSLQRNSM SPGVTHRPPS AGNTGGLMGG DLTSGAGTSA GNGYGNPRNS PGLLVSPGNL NKNMQAKSPP PMNLGMNMRK PDLRVLIPPG SKNTMPSVSQ RINNSQSAQS LATPVVSVAT PTLPGQGMGG YPSAISTTYG TEYSLSSADL SSLSGFNTAS ALHLGSVTGW QQQHLHNMP P SALSQLGACT STHLSQSSNL SLPSTQSLNI KSEPVSPPRD RTTTPSRYPQ HTRHEAGRSP VDSLSSCSSS YDGS DREDHR NEFHSPIGLT RPSPDERESP SVKRMRLSEG WAT
Source	Mammalian cell
Target Names	MEF2C
Protein Names	Recommended name: Myocyte-specific enhancer factor 2C
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.