



Recombinant Cupriavidus necator Phosphocarrier protein HPr (phbH)

Product Code	CSB-MP326411DZQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23537
Product Type	Recombinant Protein
Immunogen Species	Cupriavidus necator (strain ATCC 17699 / H16 / DSM 428 / Stanier 337) (Ralstonia eutropha)
Purity	>85% (SDS-PAGE)
Sequence	MLQRDTTIIN KLGLHARASA KLTQLAGNFV SQVKMSRNGR QVDAKSIMGV MMLAAGIGST VTLETGDPDE QEAMDALLAL IANRFGEGE
Source	Mammalian cell
Target Names	phbH
Protein Names	Recommended name: Phosphocarrier protein HPr EC= 2.7.11.- Alternative name(s): Histidine-containing protein Short name= Protein H
Expression Region	1-89
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.