



Recombinant Drosophila melanogaster Cell division control protein 2 cognate (cdc2c)

Product Code	CSB-MP326416DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23573
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MTTILDNFQR AEKIGEGTYG IVYKARSNST GQDVALKKIR LEGETEGVPS TAIREISLLK NLKHPNVVQL FDVVISGNNL YMIFEYLNMD LKKLMDKKKD VFTPQLIKSY MHQILDAVGF CHTNRILHRD LKPQNLLVDT AGKIKLADFG LARAFNVPMR AYTBEVVTW YRAPEILLGT KFYSTGVDIW SLGCIFSEMI MRRSLFPGDS EIDQLYRIFR TLSTPDETNW PGVTQLPDFK TKFPRWEGTN MPQPITEHEA HELIMSMCY DPNLRISAKD ALQHAYFRNV QHVDHVALPV DPNAGSASRL TRLV
Source	Mammalian cell
Target Names	Cdk2
Protein Names	Recommended name: Cell division control protein 2 cognate EC= 2.7.11.22 EC= 2.7.11.23
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.