



Recombinant Salmonella typhimurium DNA replication and repair protein recF (recF)

Product Code	CSB-EP326504SXB-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24900
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	SLTRLLIKD FRNIENADLA LSPGFNFLVG ANSGKTSVL EAIYTLGHGR AFRSLQPGRV IRHEQEAFVL HGRLQSEERE TSIGLTKDKQ GDSKVRIDGT DGHKIAELAH LMPMQLITPE GFTLLNGGPK YRRAFLDWGC FHNEAGFFTA WSNLKRLKQ RNAALRQVSR YEQLRPWDKE LIPLAEQIST WRAEYSSAIA QDMADTCQQF LPEFSLTFSF QRGWEKETDY ADVLERSFER DRMLTYTAHG PHKADFRIRA DGAPVEDTLS RGQLKLLMCA LRLAQGEFLT RESGRRCLYL IDDFASELDD ARRGLLASRL KATQSQVFVS AISAEHVIDM SDENSKMFTV EKGKITD
Source	E.coli
Target Names	recF
Protein Names	Recommended name: DNA replication and repair protein recF
Expression Region	2-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.