



Recombinant Streptomyces globisporus Lysozyme M1 (acm)

Product Code	CSB-MP326536SOH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25310
Product Type	Recombinant Protein
Immunogen Species	Streptomyces globisporus
Purity	>85% (SDS-PAGE)
Sequence	DTS GVQGIDVSHW QGSINWSSVK SAGMSFAYIK ATEGTNYKDD RFSANYTNAY NAGIIRGAYH FARPNASSGT AQADYFASNG GGWSRDNRTL PGVLDIEHNP SGAMCYGLST TQMRTWINDF HARYKARTTR DVVIYTTASW WNTCTGSWNG MAAKSPFWVA HWGVSAPTVP SGFPTWTFWQ YSATGRVGGV SGDVDNRNKFN GSAARLLALA NNTA
Source	Mammalian cell
Target Names	acm
Protein Names	Recommended name: Lysozyme M1 EC= 3.2.1.17Alternative name(s): 1,4-beta-N-acetylmuramidase M1
Expression Region	78-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.