



Recombinant *Saccharomyces cerevisiae* Cell cycle arrest protein BUB3 (BUB3)

Product Code	CSB-MP326630SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26449
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MQIVQIEQAP KDYISDIKII PSKSLLLITS WDGSLTVYKF DIQAKNVDLL QSLRYKHPLL CCFIDNTDL QIYVGTVQGE ILKVDLIGSP SFQALTNEA NLGICRICKY GDDKLIASW DGLIEVIDPR NYGDGVIAVK NLNSNNTKVK NKIFTMDTNS SRLIVGMNNS QVQWFRLPLC EDDNGTIEES GLKYQIRDVA LLPKEQEGYA CSSIDGRVAV EFFDDQGDDY NSSKRFAFRC HRLNLKDTNL AYPVNSIEFS PRHKFLYTAG SDGIISCWNL QTRKKIKNFA KFNEDSVVKI ACSDNILCLA TSDDTFKTNA AIDQTIELNA SSIYIIFDYE N
Source	Mammalian cell
Target Names	BUB3
Protein Names	Recommended name: Cell cycle arrest protein BUB3
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.