



Recombinant Human Vacuolar protein sorting-associated protein 37C (VPS37C)

Product Code	CSB-MP025908HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5D8V6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	METLKDKTLQ ELEELQNDSE AIDQLALESP EVQDLQLERE MALATNRS LA ERNLEFQGGL EISRSNLSDR YQELRKLVER CQEQAQKLEK FSSALQPGTL LDLLQVEGMK IEESEAMA EKFLEGEVPLE TFLENFSSMR MSLHRRVRV EKLQEVVRKP RASQELAGDA PPPRPPPPV PVPQGTTPVV EEQPQPPLAM PPYPLPYSPS PSLPVGPTAH GALPPAPFPV VSQPSFYSGP LGPTYAAQL GPRGAAGYSW SPQRSMPPRP GYPGTPMGAS GPGYPLRGGR APSPGYPQQS PYPATGGKPP YPIQPQLPSF PGQPQPSVPL QPPYPPGPAP PYGFPPPPGP AWPGY
Source	Mammalian cell
Target Names	VPS37C
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 37C Short name= hVps37C Alternative name(s): ESCRT-I complex subunit VPS37C
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.