



Recombinant Escherichia coli H repeat-associated protein yhhI (yhhI)

Product Code	CSB-MP326838ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28912
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MELKKLMEHI SIIPDYRQTW KVEHKLSLIL LLTICAVISG AEGWEDIEDF GETHLDFLKQ YGDFENGIPV HDTIARVVSC ISPAKFHECF INWMRDCHSS DDKDVAIDG KTLRHSYDKS RRRGAIHVIS AFSTMHSLVI GQIKTDEKSN EITAIPELLN MLDIKGKIIT TDAMGCQKDI AEKIQKQGGD YLFAVKGTQG RLNKAFEEKF PLKELNNPEH DSYAISEKSH GREEIRLHIV CDVPDELIDF TFEWKGLKKL CVAVSFRSII AEQKKEPEMT VRYIYSSADL TAEKFATAIR NHWHVENKLH WRLDVVMNED DCKIRRGNAE ELFSGIRHIA INILTNDKVF KAGLRRKMRK AAMDRNYLAS VLAGSGLS
Source	Mammalian cell
Target Names	yhhI
Protein Names	Recommended name: H repeat-associated protein yhhI Alternative name(s): ORF-H
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.