



Recombinant Rice gall dwarf virus Minor outer capsid protein P9

Product Code	CSB-BP326854RMJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29078
Product Type	Recombinant Protein
Immunogen Species	Rice gall dwarf virus (RGDV)
Purity	≥85% (SDS-PAGE)
Sequence	MAGKLQDGVA IAKIKETINL FCEYSFGDLV NNRREIVGRV HDARKNAALA WPDLMNCFL HSASHYGVVK FLLDIALSTR FGDFTLLGVS SQNYPFYDLH VVMTKAFCNL DFAKDEYLM I NDSFSSMMSA FLDEEGVHSA MSMELGIHDI EDRFVLRTRK LFYIIHEYHM SLDEIEPWLE KLPDASGGTL LNQKSKEQMR VIFSNAKVRI ANSINLYVTT HTNSYNEYVR EVAEYVADLW NIQTTTNTQG HENELAAEDF GVLASSSQMN GTKSELGDSV IKSDGNEVKL EPAMFTRNDD EEELAGSEFT SLLSDDGRMG
Source	Baculovirus
Protein Names	Recommended name: Minor outer capsid protein P9 Alternative name(s): Protein Pns10
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.