



Recombinant Bradyrhizobium japonicum Nitrogen fixation regulation protein fixK (fixK)

Product Code	CSB-MP326877BVW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29286
Product Type	Recombinant Protein
Immunogen Species	Bradyrhizobium diazoefficiens (strain JCM 10833 / IAM 13628 / NBRC 14792 / USDA 110)
Purity	>85% (SDS-PAGE)
Sequence	MKPSVVMIEP NGHFCSDCAI RTSAVCSSLD AAELREFEHL GRRVHFSSGE TVFSEEDITT SFYNVLEGVM RLYKLLPDGR RQIVGFALPG DFLGMNLSGR HNFSADAIGA VTVCQFAKAP FGRFIEERPQ LLRRINELAI RELSQARDHM VLLGRRSADE KVA AFL LGWR ERL LAL KGAS DTVPLPMSRQ DIADYLGLTI ETVSRTFTKL ERHGAI AIIH GGISLLDPAR VEALAAA
Source	Mammalian cell
Target Names	fixK
Protein Names	Recommended name: Nitrogen fixation regulation protein fixK
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.