



Recombinant Solanum tuberosum Actin-75 (AC75)

Product Code	CSB-BP326949FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30169
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MADVDDIQPL VCDNGTGMVK AGFAGDDSPR AVFPSIVGRP RHTGVMVGMG QKDAYVGDEA QSKRGILTWK YPIEHGIVNN WDDMEKIWHH TFYNELRVAP EEHPVLLTEA PLNPKANREK MTQIMFETFN XXXXXXXXXXXX XXXXXXXXXX XXSVLDSGDG VSHTVPIYEG YALPHAILRL DLAGRDLTEH LAKILTERGY SFTTSAEKEI VRDMKEKLAY IALDFEQELD MVKSSSSVEK NFELPDGQVI TIGAERFRCP EVLFQPSTIG MEAAGIHETT YNSIMKCDVD IRKDLYGNIV LSGGSTMFPG IADRMSKEIT ALAPSSMKIK VVAPPERKYS VWIGGSILAS LSTFQQMWIA KAEFDESGPS IVHRKCF
Source	Baculovirus
Target Names	AC75
Protein Names	Recommended name: Actin-75
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.