



Recombinant Human Coiled-coil domain-containing protein 151 (CCDC151)

Product Code	CSB-YP004670HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5D8V7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MTSPLCRAAS ANALPPQDQA STPSSRVKGR EASGKPSHLR GKGTAQAWTP GRSKGGSFHR GAGKPSVHSQ VAELHKKIQL LEGDRKAFFE SSQWNIKKNQ ETISQLRKET KALELKLLDL LKGDEKVVQA VIREWKWEKP YLKNRTGQAL EHLDHRLREK VKQQNALRHQ VVLRQRRLEE LQLQHSLRLL EMAEAQNRHT EVAKTMRNLE NRLEKAQMKA QEAHITSVY LQLKAYLMDE SLNLENRLDS MEAEVVRTKH ELEALHVVNQ EALNARDIAK NQLQYLEETL VRERKKRERY ISECKKRAEE KKLNERMER KTHREHLLLQ SDDTIQDSLH AKEEELRQRW SMYQMEVIFG KVKDATGTDE THSLVRRFLA QGDTFQALET LKSENEQTLV RLKQEKQQLQ RELEDLKYSG EATLVSQQKL QAEAQERLKK EERRHAEAKD QLERALRAMQ VAKDSLEHLA SKLIHITVED GRFAGKELDP QADNYVPNLL GLVEEKLLKL QAQLQGHDVQ EMLCHIANRE FLASLEGRLP EYNTRIALPL ATSKDKFFDE ESEEDNEVV TRASLKIRSQ KLIESHKKHR RSRRS</p>
Source	Yeast
Target Names	CCDC151
Protein Names	Recommended name: Coiled-coil domain-containing protein 151
Expression Region	1-595
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.