



Recombinant Escherichia coli Putative prophage Qin DNA-packaging protein NU1 homolog (nohA)

Product Code	CSB-BP327019ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31061
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MEVNKKQLAD IFGASIRTIQ NWQEQGMPVL RGGGKGNEVL YDAAVIRWY AERDAEIE NE KLRREVEELR QASETDLQPG TIEYERHRLT RAQADAQELK NARDSAEVVE TAFCTFVLSR IAGEIASILD GIPLSVQRRF PELENRHVDF LKRDIKAMN KAAALDELIP GLLSEYNRAD RQYAGGSR S
Source	Baculovirus
Target Names	nohA
Protein Names	Recommended name: Putative prophage Qin DNA-packaging protein NU1 homolog
Expression Region	1-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.