



Recombinant Human Proteasome subunit beta type-11 (PSMB11)

Product Code	CSB-BP018878HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5LHX3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	T TTLAFRFRHG VIAADTRSS CGSYVACPAS CKVIPVHQHL LGTTSGTSAD CATWYRVLQR ELRLRELREG QLPSVASAAK LLSAMMSQYR GLDLCVATAL CGWDRSGPEL FYVYSDGTRL QGDIFSVGSG SPYAYGVLDLDR GYRYDMSTQE AYALARCAVA HATHRDAYS GSVDLFHVRE SGWEHVSRS D ACVLYVELQK LLEPEPEEDA SHAHPEPATA HRAAEDRELS VGPGEVTPGD SRMPAGTETV
Source	Baculovirus
Target Names	PSMB11
Protein Names	Recommended name: Proteasome subunit beta type-11 EC= 3.4.25.1 Alternative name(s): Proteasome subunit beta-5t
Expression Region	50-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.