



# Recombinant *Saccharomyces cerevisiae* Protein MLF3 (MLF3)

<b>Product Code</b>	CSB-EP327103SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P32047
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MCVYKSNSNN SNPSFIFERT VQEASSNDLF LQPPVSASNT SHSSRSNSFY NLQTISPIPI SGSEVRTPSL RKNSNNVSSP LDNVIPTSRAS NSTTSSLA HQEYILNPIC NMQNHHRRR TLENSVAPAL DASCIVNDE NTDLSVDVMV YSRRPSSAVS LNMALLARTN SATLPSSSESS PASPDLKLSR SHSHAATRP TLNNINNTGM TTTTNGEPN SRILRFYSYV DMLNDEKLAQ ANNTPTSRPP MKSQAYSCPF ILKRSPQAY SSSSATTTFS NPFIKTTELP ATSPYVSPQQ SARQYSNNAN NNAKSPKNRS SSILFQRQSI LSNVDPVANM HKNPKFQIES SDSEEDLTM DMLDPSFPLS SSLRSSANLA SNP ELATQTP LSTSSSYTAI GKPMPLSTDP SYVSSSNTLS SEHEL RVEKV SEVLKKKVS N GGFSTEFNSC DT
<b>Source</b>	E.coli
<b>Target Names</b>	MLF3
<b>Protein Names</b>	Recommended name: Protein MLF3 Alternative name(s): Multicopy suppressor of leflunomide sensitivity 3
<b>Expression Region</b>	1-452
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.