



Recombinant Enterobacteria phage T4 RNA ligase 2 (Y10A)

Product Code	CSB-BP327123EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32277
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MFKKYSSLEN HYN SKFIEKL YSLGLTGGEW VAREKIHGTN FSLIIERDKV TCAKRTGPIL PAEDFFGYEI ILKNYADSIK AVQDIMETSA VVS YQVFGEF AGPGIQKNVD YCDKDFYVFD IIVTTESGDV TYVDDYMMES FCNTFKFKMA PLLGRGKFEE LIKLPNDLDS VVQDYNFTVD HAGLVDANKC VVNAEAKGEV FTAEGYVLKP CYP SWLRNGN RVAIKCKNSK FSEKKKSDKP IKAKVELSEA DNKLVGILAC YVTLNRVNNV ISKIGEIGPK DFGKVMGLTV QDILEETSRE GITLTQADNP SLIKKELVKM VQDVL RP AWI ELVS
Source	Baculovirus
Target Names	Y10A
Protein Names	Recommended name: RNA ligase 2 Short name= Rnl2 EC= 6.5.1.3
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.