



Recombinant *Saccharomyces cerevisiae* Histone chaperone ASF1 (ASF1)

Product Code	CSB-EP327134SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32447
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSIVSLLGIK VLNNPAKFTD PYEFEITFEC LESLKHDLW KLTYVGSSRS LDHDQELDSI LVGPVPVGVN KVFVSADPPS AELIPASELV SVTVILLSCS YDGRFVVRVG YYVNNEYDEE ELRENPPAKV QVDHIVRNIL AEKPRVTRFN IVWDNENEGD LYPPEQPGVD DDEEEDDEEE DDEDEDDEDE DDDQEDGEGE AEEAAEEEEEE EEEKTEDNET NLEEEEEEDIE NSDGDEEEGE EEVGSVDKNE DGNDKKRRKI EGGSTDIEST PKDAARSTN
Source	E.coli
Target Names	ASF1
Protein Names	Recommended name: Histone chaperone ASF1 Alternative name(s): Anti-silencing function protein 1 Short name= yASF1
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.