



Recombinant *Saccharomyces cerevisiae* Trans-aconitate 3-methyltransferase (TMT1)

Product Code	CSB-EP327154SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32643
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	STFSASDFN SERYSSSRPS YPSDFYK MID EYHDGERKLL VDVGCGPGTA TLQMAQELKP FEQIIGSDL S ATMIKTAEVI KEGSPD TYKN VSFKISSDD FKFLGADSVD KQKIDMITAV ECAHWDFEK FQRSAYANLR KDGTIAIWGY ADPIFPDYPE FDDLMI E VPY GKQGLGPYWE QPGRSRLRNM LKDSHLDPEL FHDIQVSYFC AEDVRDKVKL HQHTKKPLLI RKQVTLVEFA DYVRTWSAYH QWKQDPKNKD KEDVADWFIK ESLRRRPELS TNTKIEVVWN TFYKLGKRV
Source	E.coli
Target Names	TMT1
Protein Names	Recommended name: Trans-aconitate 3-methyltransferase EC= 2.1.1.145
Expression Region	2-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.