



Recombinant Burkholderia cepacia Phthalate dioxygenase reductase (ophA1)

Product Code	CSB-YP327204BXR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33164
Product Type	Recombinant Protein
Immunogen Species	Burkholderia cepacia (Pseudomonas cepacia)
Purity	>85% (SDS-PAGE)
Sequence	TTPQEDGFL RLKIASKEKI ARDIWSFELT DPQGAPLPPF EAGANLTVAV PNGSRRTYSL CNDSQERNRY VIAVKRDSNG RGGISFIDD TSEGDAVEVS LPRNEFPLDK RAKSFILVAG GIGITPMLSM ARQLRAEGLR SFRLYYLTRD PEGTAFFDEL TSDEWRSDVK IHHDHGDPTK ADFDWSVFEK SKPAQHVYCC GPQALMDTVR DMTGHWPSGT VHFESFGATN TNARENTPFT VRLSRSGTSF EIPANRSILE VLRDANVRVP SSCESGTCGS CKTALCSGEA DHRDMVLRDD EKGTQIMVCV SRAKSAELVL DL
Source	Yeast
Target Names	ophA1
Protein Names	Recommended name: Phthalate dioxygenase reductase Short name= PDR EC= 1.-.-
Expression Region	2-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.