



Recombinant *Drosophila sechellia* Accessory gland-specific peptide 26Aa (Acp26Aa)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-BP327253DMG |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P33736 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila sechellia</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | NC DGEHQLDSSM DLKSDSTKSA VLKNVAPKND ATQAEIAKDD VALKSGKKGD YVMEIDVSDI PLDDYPINNS KSRKNSSTLP SQILTDKPNQ GSNQIALKAL KHRLLEMEQNN NLFLRNHSVS LMNEIEARKT DIIQARQLNI DLELELEALK RELSEMNAQN ARKSTKSCKK RPSKDIASPV NQLQEIVVKN TYRNKYL TLL TQLAQKINYE IANVNNPATD VPTGKSPSEG NPSTT |
| Source | Baculovirus |
| Target Names | Acp26Aa |
| Protein Names | Recommended name: Accessory gland-specific peptide 26Aa Alternative name(s): Male accessory gland secretory protein 355A |
| Expression Region | 19-255 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |