



Recombinant Drosophila sechellia Accessory gland-specific peptide 26Aa (Acp26Aa)

Product Code	CSB-YP327253DMG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33736
Product Type	Recombinant Protein
Immunogen Species	Drosophila sechellia (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	NC DGEHQLDSSM DLKSDSTKSA VLKNVAPKND ATQAEIAKDD VALKSGKKGD YVMEIDVSDI PLDDYPINNS KSRKNSSTLP SQILTDKPNQ GSNQIALKAL KHROLLMEQNN NLFLRNHSVS LMNEIEARKT DIIQARQLNI DLELELEALK RELSEMNAQN ARKSTKSCKK RPSKDIASPV NQLQEIVIVKN TYRNKYLLTLL TQLAQKINYE IANVNNPATD VPTGKSPSEG NPSTT
Source	Yeast
Target Names	Acp26Aa
Protein Names	Recommended name: Accessory gland-specific peptide 26Aa Alternative name(s): Male accessory gland secretory protein 355A
Expression Region	19-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.