



Recombinant Escherichia coli Uncharacterized protein yejH (yejH)

Product Code	CSB-EP327272ENV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33919
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MIFTLRPYQQ EAVDATLNHF RRHKTPAVIV LPTGAGKSLV IAELARLARG RVLVLAHVKE LVAQNHAKYQ ALGLEADIFA AGLKRKESHG KVVFGSVQSV ARNLDAFQGE FSLIVDECH RIGDDEESQY QQILTHLTKV NPHLRLGLT ATPFRLGKGW IYQFHYHGMV RGDEKALFRD CIYELPLRYM IKHGylTPPE RLDMPVVQYD FSRLQAQSNQ LFSEADLNRE LKKQQRITPH IISQIMEFAA TRKGVMIFAA TVEHAKEIVG LLPEDAALI TGDTPGAERD VLIENFKAQR FRYLVNVAVL TTGFDAPHVD LIALRPTES VSlyQQIVGR GLRLAPGKTD CLILDYAGNP HDLYAPEVGT PKGKSDNVPV QVFCPACGFA NTFWGKTTAD GTLIEHFGRR CQGWFEEDDG HREQCDFRFR FKNCPQCNAE NDIAARRCRE CDTVLVDPDD MLKAALRLKD ALVLRCSGMS LQHGHDEKGE WLKITYYDED GADVSERFRL QTPAQRtAFE QLFIRPHTRT PGIPLRWITA ADILAQQALL RHPDFVARM KGQYWQVREK VFDYEGRFRL AHELRG
Source	E.coli
Target Names	radD
Protein Names	Recommended name: Uncharacterized protein yejH
Expression Region	1-586
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.