



Recombinant Uncharacterized protein B0303.3 (B0303.3)

Product Code	CSB-MP327303CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34255
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MLRAVSTSFG TARAASAVAK KNMPNIVLVD AVRTPFVVSF TVFKDLMAVD LQKEAIKALV ETKLTPYEQL DHICGTVIQ ECKTSNIARE AALLAGVPDK IPAHTVTLAC ISSNVAMTTG MGMLATGNAN AIIAGGVELL SDVPIRYNRN ARKAMLG MNK AKDVPSKLKI GGQIVKNLLS PELPAVAEFS TGETMGHSGD RLAAAFNVS R REQDEFAIRS HTLASEAAKN GKFTDVVPVF LDGKKPKTIK EDNGIRVSTL EKLSSLKPAF VKPHGTVTAA NASYLT DGAS AALIMTEEYA LANGYKPKAY LRDYLYVAQD PKDQLLLSPA YVIPKLLDKA GLTLKDVDVF EIHEAFAGQV LANLNAMDS YFCKEQMKRS GKFG RVPMDK LNLWGGSLSI GHPFGATGVR LATHSAHRLK EEKGQYAVIA ACAAGGHGVG MLIEAYGK
Source	Mammalian cell
Target Names	B0303.3
Protein Names	Recommended name: Uncharacterized protein B0303.3
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.