



Recombinant Uncharacterized protein F54H12.2 (F54H12.2)

Product Code	CSB-YP327328CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34456
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPGTQTCIQN SRWTVVNLKN AFQAEGPWEF VLGNNRSYLNKRTYLVYTFNITDQNGAI VNIPNPITKD SLVYAPINNI AHSIVKNFSL HINSQLAYHNSSNYAYKAYF ENALMYGKEI KNSTLTAAGY FHEDEIGTPT SKGFLNRCGTGTTQVAANIS LDLMNQPRVL LNSCNVKLTA YPNNSEFLIE AYNHGQQKFKFNVTDVYALV NEFDLTDGLS NELEAAIEH KMIQYPMISS QVRSFYIDGG RFDAPANTLF TSKMPRRLFL GLVSSEAYNG SFGTTPFNFK PYGVTDVHIDYCGQTIPGRP MNLDFENKVF VEAYVQLQET LGHTRNNFSC NSIDIGMFRSKGYTIFGFEL SPVAQDNNLF ELVRQTNVSV RLNFKKETPP GGLYCVVYAEFDQIFSLDFM RNPIVDSIV
Source	Yeast
Target Names	F54H12.2
Protein Names	Recommended name: Uncharacterized protein F54H12.2
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.