



Recombinant Zinc finger-interacting factor 1 (zif-1)

Product Code	CSB-YP327331CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34482
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSECSASTSQ LSTFCVNEIP WKLVPQTCDE FSEYEFIDEK YRRLSDVESF VWKWSPCGRY GIDLPGGRDV INQQLRDRMA VFDLHTNKWT LYTMDRGAFE IPRVCFLYWA RNDIIGTITL VKDYEAGPDD EIPLKFKQSF FYVSHETQTL VYSGSCKYRT NKTIVASYWY HIFENRRGEP FIDHNNYLWL VFSNGVDSL ILLPFIPTSTC FKA KCVSFDM NDI VSQFTGI RVFLDPFVSS HYVPWVYRND V NFFIRIDSH SFIEEQNRDY GLDPQIDNPE TKPHRGVFQL RVNLKDVI EK GEKDAIHAVN PETIETRFNR GDIELWRKIS DELRQEHQSI AQHGKMMVVVQ RWSHFTVDMG LSKKMRTRVA DLDARDYCCS HKPRNSSFCC IDLDSERIKA VDALILGTDV VLPHPSGSVF LFRYKENYGM SYCELPFFQP LSLKLCIQV LKESVEIDPS YYRIKCCCLKM QEELEKETAE NLCTINIQQ
Source	Yeast
Target Names	zif-1
Protein Names	Recommended name: Zinc finger-interacting factor 1 Short name= Protein zif-1
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.