



Recombinant Probable glutamine synthetase (gln-2)

Product Code	CSB-MP327333CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34497
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTHLNFETRM PLGTAVIDQF LGLRPHPTKI QATYVWIDGT GENLRSKTRT FDKLPKRIED YPIWNYDGSS TGQAKGRDSD RYLRPVAAYP DPFLGGGNKL VMCDTLDHQM QPTATSHRQA CAEIMHEIRD TRPWFGEQEQE YLIVDRDEHP LGWPKHGFPD PQGKYYCSVG ADRAFGREVV ETHYRACLHA GLNIFGTNAE VTPGQWEFQI GTCEGIDMGD QLWMSRYILH RVAEQFGVCV SLDPKPKVTM GDWNGAGCHT NFSTAEMRAP GGIAAIEAAM TGLKRTHLEA MKVYDPHGGE DNLRRLTGRH ETSSADKFSW GVANRGCSIR IPRQVAAERK GYLEDRRPSS NCDPYQVTAM IAQSILF
Source	Mammalian cell
Target Names	gln-2
Protein Names	Recommended name: Probable glutamine synthetase EC= 6.3.1.2 Alternative name(s): Glutamate--ammonia ligase
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.