



Recombinant Uncharacterized protein R10E12.2 (R10E12.2)

Product Code	CSB-BP327340CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34553
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MRGALLLLVL LSAQHVSTSQ FWLSVEVEQI DWTGCLTTA LCSHPRFQLI KDLLPVNEKV TMNWPIVEHF DKESHRPFVS YWPSGRTEI SMSAQVVGTD RTYGFPRICD QSPSVRIFPE EHKKIVAHLE KEKPTGKPPM DSLKIKVKGK CFNATMTVIK HTERCPWCPD PKEITIIGQE PGSEATGLRA GSAAWLFGSS SSLISDDRIV HIGVLVLAIV AVLSSTAFAI ILVMYLRNKR LVKETLKKPR FHPYISVKGH EIAEDNCRYD LPWEQQQRPL TYWMTSSNKS SESTMTSPLD SASSLHGSGG SHQNPHMYNS HHHQHQQQQN AHPSEMYHST YTRDGYQTYR PPPPSVHPPI FPPQTQLFHP PTYSTQRHVT SPNSSRNDDS GLESV
Source	Baculovirus
Target Names	R10E12.2
Protein Names	Recommended name: Uncharacterized protein R10E12.2
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.