



Recombinant Uncharacterized cyclin-like protein ZK353.1 (ZK353.1)

Product Code	CSB-EP327349CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34624
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MGNSSCCLRT RSSSGEDKSY NNDGQYIRTN QVEFQYVNQV FPRDETSTNF LPHISEREVT EGYEEDPSTN PTARPTFMER SKSEMCLKDN RRSCYMLDAL AAGGHHPGIL PRSLRKSSSC STIYIDDSTV SQPHLKNTIK CISLAIYYHI SNRKNRGHER LMEIFEERLH PIFRDPIPE QMTRDPDHRN IYRFVRNLFS SAQLTAECAI ITLVYIERLL NYAEMDLCPN NWRRVVLGSI MLASKVWDDQ AVWNVDYCQI LRDTNVDDMN ELERRFLECL DFNIEVPSSV YAKYYFDLRT LALANDLQLP IQPLYKERAQ RLEALSRVFE DKIQSSSLPK RARSAEHLVF EHPAVLS
Source	E.coli
Target Names	cyt-1
Protein Names	Recommended name: Uncharacterized cyclin-like protein ZK353.1
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.