



Recombinant Leucine aminopeptidase 1 (lap-1)

Product Code	CSB-MP327350CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34629
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTQVLVRNGI QAVGDGLTSL IIVGKKSVLK NVTFEGKFKE VAQKFVTDGD SWNSMISRIP ASGRHPLHYE LAHLITVPDA SSRGNTPTNA HSIYKELKPI NYPEDTKNVH FVLFAEYPDV LSHVAAIART FCKFSMKTSG IRELVNIDV VCDKLTNEDA VFLTDLSESV RETARLIDTP ANILTTDALV DEAVKVG NAT GSKITVIRGE ELLKAGFGGI YHVGKAGPTP PAFVVL SHEV PGSTEHIALV GKGVVYDTGG LQIKTKTGMP NMKRDMGGAA GMLEAYSALV KHGFSQTLHA CLCIVENNVS PIANKPDDII KMLSGKTVEI NNTDAEGR LI LADGVFYAKE TLKATTIFDM ATLTGAQAWL SGR LHGAAMT NDEQLENEII KAGKASGDLV APMLFAPDLF FGDLKSSIAD MKNSNLGKMD GPPSAVAGLF IGAHIGFGEG LRWLHLDIAA PAEVGDRGTG YGPALFSTLL GKYTSVPMLK Q
Source	Mammalian cell
Target Names	lap-1
Protein Names	Recommended name: Leucine aminopeptidase 1 EC= 3.4.11.1
Expression Region	1-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.