



Recombinant *Saccharomyces cerevisiae* Protein UTH1 (UTH1)

Product Code	CSB-YP327496SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36135
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	APA VHSDNHHHN DKRAVVTVTQ YVNADGAVVI PAATTATSAA ADGKVESVAA ATTTLSSTAA AATTSAAASS SSSSSSSSSS SSSVGSGDFE DGTISCSDFP SGQGAVSLDW LGLGGWASIM DMNGNTATSC QDGYCYAC SPGYAKTQWP SEQPSDGRSV GGLYCKNGKL YRSNTDTNSL CVEGQGSAQA VNKVSGSIAI CGTDYPGSEN MVVPTVVGAG SSQPINVIKE DSYYQWQGKK TSAQYYVNNA GVSVEDGCIW GTEGSGVGNW APVVLGAGYT DGITYLSIIP NPNNKEAPNF NIKIVATDGS TVNGACSYEN GVYSGSGSDG CTVSVTSGSA NFVIFY
Source	Yeast
Target Names	UTH1
Protein Names	Recommended name: Protein UTH1 Alternative name(s): Youth protein 1
Expression Region	18-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.