



# Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YKR075C (YKR075C)

<b>Product Code</b>	CSB-EP327498SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P36155
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTSLDDTIIS YQNIMLLDNM TNYNKPAIDY FHHEFNDA SL EISASWTLLL KMRKHKLLRL PSCSSEVDLD YNMYLVRLHH CLWRRWSINH YGLQNSKSNP LSINWVKETD VTVLYGPDLT NIDSNENEIS PVQIQIDQKQ TKNLKSALKK NTECWVTEEV DEINASIESN DNALVKLEDI SCPSSVDSHT SSIFDQHSTC TKISSIDEDS EDLMNEKKEQ FPRKLKFNQA VMKREIDSKG TIRESLININ DIQHSRHHRR HRRRHHHHHH QNSSHSDETI KEAHYEF SNY TFGTMEEDIF YRNQVVF
<b>Source</b>	E.coli
<b>Target Names</b>	YKR075C
<b>Protein Names</b>	Recommended name: Uncharacterized protein YKR075C
<b>Expression Region</b>	1-307
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.