



Recombinant Schizosaccharomyces pombe Fructose-bisphosphate aldolase (fba1)

Product Code	CSB-BP327559SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36580
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MGILDIVPTG VITGDNVLKL FTYAREHGFA IPAINVTSSS TAIAALEAAR EARSPIILQT SNGGAHFFAG KESSNEGQKA SIAGSIAAAH YIRSIAPFFG VPVVMHSDHC AKKLLPWMDG MFEADEAYFK IHGEP LFSSH MLDLSEEPKK ENIAQVKEYC KRAVPMKIWI EMEIGITGGE EDGVDNSHVS HTELYTQPED IWDVYRELSS VTPYFSIAAA FGNVHG VYKP GNVK LQPALL GQH QAYVKEQ LKT TNDKPVF FVFHGGSGSS VNEFRTGIKC GVVKNIDTD TQFAYVEGVR DYVLKYKDYL MTPVGNPEGA DKPNKKKFDP RVWIHEGEKT MTKRVLTALE DFYTVNTL
Source	Baculovirus
Target Names	fba1
Protein Names	Recommended name: Fructose-bisphosphate aldolase Short name= FBP aldolase Short name= FBPA EC= 4.1.2.13 Alternative name(s): Fructose-1,6- bisphosphate aldolase
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.