



Recombinant Schizosaccharomyces pombe Serine/threonine-protein kinase csk1 (csk1)

Product Code	CSB-BP327563SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36615
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKSVGHFVPW LDIRHLTDG TISEVFGGER KNSKKLYVIK VQGLVFKRPP HDAMRGKLIL ESIGHPIER IVDSFIDNEA GSVYLITFSK SFVLSVMDDE ISIDTKCKIV LQISSALEYL EKHGILHRDI HPNNILLDSM NGPAYLSDFS IAWSKQHPGE EVQELIPQIG TGHYRAIETL FGCHSYGHEV DRWTFGILIA ELFSNQALFD DGSSEGWPSE LRLTSSIIQT LGTPNPSMWP ELSTFPDWNK FIFHEYPPKP WSEILPSVDT SIQYIVSHLV TYSNRASPSF VIESFPKVSA RLSQYA
Source	Baculovirus
Target Names	csk1
Protein Names	Recommended name: Serine/threonine-protein kinase csk1 EC= 2.7.11.22 Alternative name(s): CAK-activating kinase Short name= CAKAK
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.