



# Recombinant *Saccharomyces cerevisiae* Low specificity L-threonine aldolase (GLY1)

<b>Product Code</b>	CSB-BP327657SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P37303
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTEFELPPKY ITAANDLRSD TFFTPTAEMM EAALEASIGD AVYGEDVDTV RLEQTVARMA GKEAGLFCVS GTLSNQIAIR THLMQPPYSI LCDYRAHVYT HEAAGLAIS QAMVVPVPS NGDYLTLEDI KSHYVPDDGD IHGAPTRLIS LENTLHGIVY PLEELVRIKA WCMENGLKLH CDGARIWNAA AQSGVPLKQY GEIFDSISIC LSKSMGAPIG SVLVGNLKFV KKATHFRKQQ GGGIRQSGMM ARMALVNINN DWKSQLLYSH SLAHELAEYC EAKGIPLESP ADTNFVFINL KAARMDPDVL VKKGLKYNVK LMGGRVSFHY QVTRDTLEKV KLAISEAFDY AKEHPFDCNG PTQIYRSEST EVDVDGN AIR EIKTYKY
<b>Source</b>	Baculovirus
<b>Target Names</b>	GLY1
<b>Protein Names</b>	Recommended name: Low specificity L-threonine aldolase Short name= Low specificity L-TA Short name= TA EC= 4.1.2.48
<b>Expression Region</b>	1-387
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.