



Recombinant Escherichia coli Uncharacterized HTH-type transcriptional regulator yiaU (yiaU)

Product Code	CSB-YP327714ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37682
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MTKLQLKYRE LKIISVIAAS ENISHAATVL GIAQANVSKY LADFESKVGL KVFDRTTRQL MLTPFGTALL PYINDMLDRN EQLNFIADY KHEKRGRVTI YAPTGIITYL SKHVIDKIKD IGDITLSLKT CNLERNAFYE GVEFPDDCDV LISYAPPKDE SLVASFITQY AVTAYASQRY LEKHPISRPD ELEHHSCILI DSMMIDDANI WRFNVAGSKE VRDYRVKGNV VCDNTQSALE LARNHLGIVF APDKSVQSDL QDGTLVPCFQ QPYEWWLDLV AIFRKREYQP WRVQYVLDEM LREIRHQLAQ SQQLRPEQAA ESED
Source	Yeast
Target Names	yiaU
Protein Names	Recommended name: Uncharacterized HTH-type transcriptional regulator yiaU
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.